

ABISee, Inc.
Zoom-Office Quick Reference Guide (3.0)

Keystrokes for reading functions.

Start/Stop Zoom-Ex Camera: Space or (Caps Lock + Dot) .

Start reading: E or (Caps Lock +Down) .

Pause/Resume reading: X or (Caps Lock +End) .

Reading by character, word, sentence or section: use standard Windows commands with arrow keys.

Cancel processing: Escape.

Speech output control keystrokes.

Increase voice speed: (Shift + Q) or F9.

Decrease voice speed: (Shift + Z) or F8.

Hear page description: Shift + E.

Hear current control: Shift.

Stop speaking: Control.

ABISee, Inc.

Zoom-Office Quick Reference Guide (3.0)

Zoom-Ex Visual Keyboard Command Keys.

Zoom In/Out : + / - .

**Switch between Preview/Magnification and Formatted modes:
F2 (or F).**

Go one screen Up/Down in reformatted mode: Page Up / Down

Toggle color: F3 (or C) .

Toggle sharpness: F4 (or S).

Rotate image 90 degrees: F5 (or R) .

Show/Hide preview panel in magnified mode: F6 (or V) .

Adjust scroll speed in Row mode: . / , (or < / >) .

To simulate X-Y table use arrow keys. (Note:“Use arrow keys to simulate X-Y table option” must be enabled).

While camera is ON you can:

Turn LED On / Off: Control + L.

Adjust brightness: . / , (or < / >) .

Book Keyboard Command Keys

Add page: A (while creating a new book).

Next page: N.

Previous page: P.